

Medical & Scientific Affairs

Are 'patient-ready' endoscopes free of microbial contamination?

- Gastrointestinal endoscopies are high volume procedures performed worldwide.¹⁻⁵
- Researchers have linked inadequate reprocessing to outbreaks involving multidrug-resistant bacteria.⁶
- Bacterial contamination continues to persist despite reprocessing endoscopes according to guidelines and endoscope instructions for use.^{1,6,7}
- High-level disinfection may be unable to completely remove biofilm originating from residual bacterial contamination.⁸⁻¹²



Figure 1.
Location of studies reporting endoscope-associated infections from inception to 2020.

Adapted from Deb et al.¹

● Duodenoscopy and ERCP ■ Upper GI Endoscopy ▲ Colonoscopy and Sigmoidoscopy

Every year, endoscopies are regularly performed around the world^{1,3-5}

Annually, health care professionals perform millions of gastrointestinal (GI) endoscopy procedures around the world.^{1,3-5} An estimated 17.7 million GI endoscopies occur each year in the United States.^{1,2} In Europe, the annual number of procedures is in the tens of millions.³ In Japan and China, more than 14 million procedures are performed every year.^{4,5}






After each procedure, reprocessing staff are responsible for cleaning and disinfection of endoscopes.⁶ Reprocessing lapses result in contaminated endoscopes which put patients at risk of infections.⁹

Contaminated GI endoscopes implicated in more outbreaks than other medical devices^{6,9}

Based on a review of 15 journal articles, the proportion of duodenoscope-associated infected or colonised patients ranged from 6% to $\geq 20\%$.⁷ From infection rates reported in 16 studies, the calculated composite infection rate, which included duodenoscopes and gastroscopes, was 123 per 1,000 procedures.¹ An emerging cause of endoscope-associated infections is multidrug-resistant organisms.^{1,6,7}

Reports of more than 500 episodes of microbial transmission span 45 years in a review of 63 articles. In five outbreaks with no reprocessing breaches, bacterial transmission resulted in 93 infected patients.¹⁰ Post-procedure infections arise from contaminated endoscopes or exposure to the patient's own gut flora.¹

Complications from contaminated endoscopes^{1,7}

-  Intestinal colonisation
-  Delayed remote site infections
-  Secondary transmission to patients at other hospitals
-  Bloodstream infections, sepsis
-  Deaths

Contamination persists despite endoscopes reprocessed according to guidelines^{1,6,7}

Almost three-quarters of endoscopes sampled (8/11; 73%) were still contaminated after a successful high-level disinfection (HLD) cycle was completed. Additionally, surface ATP was detected on two endoscopes, while surface protein was present on six endoscopes.¹³

Nearly half of reprocessed GI endoscopes (47/102; 46.1%) were found contaminated in an Italian teaching hospital. These endoscopes were positive for *E.coli* (of which one was multidrug resistant), *K. pneumoniae* and multidrug-resistant *P. aeruginosa*.¹⁴ Other researchers detected biofilm contamination in endoscopes reprocessed per guidelines.¹⁵

Suggestions of biofilm found in endoscope deemed source of outbreak¹⁶


An outbreak involved three patients infected with multidrug-resistant *P. aeruginosa* sepsis. All patients underwent endoscopic retrograde cholangiopancreatography (ERCP) with the same endoscope. After intensive HLD, negative cultures suggested the endoscope was patient-ready.¹⁶

Patients and endoscope channels yielded linked *P. aeruginosa* isolates over several months. Four months after ethylene oxide sterilisation, *P. aeruginosa* contamination reoccurred. Manufacturer repair of the endoscope found suggestions of biofilm inside the endoscope channels. Persistent contamination after HLD and sterilisation highlight the difficulty of removing biofilm.¹⁶

The challenges of biofilm removal^{8-12,17,18}

The biofilm in endoscopes forms under multiple cycles of wetting and drying. This cyclic buildup of biofilm results in compacted biofilm which is difficult to remove.¹⁹ Researchers in one study observed damage, residue or debris on or inside all reprocessed endoscopes. These defects may harbour bioburden and could facilitate biofilm formation.¹⁷

Rapid biofilm formation may occur in a new endoscope after only 30 days of clinical practice.¹⁸ Limited access for brushes, water or air contributes to ineffective reprocessing.¹⁰⁻¹² Evidence shows HLD is sometimes unable to completely remove biofilm.^{9,9}





[Click here](#) to learn more about endoscope reprocessing



Nanosonics Limited (Manufacturer)

7-11 Talavera Road, Macquarie Park, NSW 2113, Australia
T: +61 2 8063 1600 E: info@nanosonics.com.au
W: www.nanosonics.com.au

Nanosonics UK Limited (Distributor)

Unit 2, Linfit Court, Colliers Way, Clayton West, Huddersfield, HD8 9WL, UK
T: +44 1484 860581 E: info@nanosonics.co.uk W: www.nanosonics.co.uk

Nanosonics Europe GmbH Ireland (Distributor)

Unit 1, Building 5, Port Tunnel Business Park
Clonsbaugh, Dublin 17, D17 P497, Ireland
T: +353 (0)1 6835790 E: irlinfo@nanosonics.com W: www.nanosonics.ie

Nanosonics Europe GmbH (EU Representative)

Poppenbütteler Bogen 66, 22399 Hamburg, Germany
T: +49 40 46856885 E: info@nanosonics.eu W: www.nanosonics.eu

Nanosonics France (Distributor)

2, route de la Noue 91190, Gif sur Yvette, France
T: +33 1 64 86 58 59 E: info@nanosonics.fr W: www.nanosonics.fr

ATP: adenosine triphosphate; ERCP: endoscopic retrograde cholangiopancreatography; GI: gastrointestinal; HLD: high-level disinfection.

References

1. Deb A, Perisetti A, Goyal H, et al. Gastrointestinal Endoscopy-Associated Infections: Update on an Emerging Issue. *Dig Dis Sci*. 2022;67(5):1718-1732.
2. Peery AF, Crockett SD, Murphy CC, et al. Burden and Cost of Gastrointestinal, Liver, and Pancreatic Diseases in the United States: Update 2018. *Gastroenterology*. 2019;156(1):254-272.e11.
3. Rutter MD, Senore C, Bisschops R, et al. The European Society of Gastrointestinal Endoscopy Quality Improvement Initiative: developing performance measures. *Endoscopy*. 2016;48:81-89.
4. Tanaka K. Japan Endoscopy Database project. *Dig Endosc*. 2022;34 Suppl 2:20-22.
5. Zhou S, Zhu Z, Dai W, et al. National survey on sedation for gastrointestinal endoscopy in 2758 Chinese hospitals. *Br J Anaesth*. 2021;127:56-64.
6. Grein JD, Murthy RK. New Developments in the Prevention of Gastrointestinal Scope-Related Infections. *Infect Dis Clin N Am*. 2018;32(4):899-913.
7. Ofstead C, Buro B, Hopkins K, Eiland J, Wetzler H, Lichtenstein D. Duodenoscope-associated infection prevention: A call for evidence-based decision making. *Endosc Int Open*. 2020;08(12):E1769-E1781.
8. Brunke MS, Konrat K, Schaudinn C, et al. Tolerance of biofilm of a carbapenem-resistant *Klebsiella pneumoniae* involved in a duodenoscopy-associated outbreak to the disinfectant used in reprocessing. *Antimicrob Resist Infect Control*. 2022;11(1):81.
9. Lichtenstein D, Alfa MJ. Cleaning and Disinfecting Gastrointestinal Endoscopy Equipment. *Clin Gastrointest Endosc*. 2019;(JAMA 312 2014):32-50.e5.
10. Kovaleva J. Infectious complications in gastrointestinal endoscopy and their prevention. *Best Pract Res Clin Gastroenterology*. 2016;30(5):689-704.
11. Rauwers AW, Troelstra A, Fluit AC, et al. Independent root-cause analysis of contributing factors, including dismantling of 2 duodenoscopes, to investigate an outbreak of multidrug-resistant *Klebsiella pneumoniae*. *Gastrointest Endosc*. 2019;90(5):793-804.
12. Ribeiro MM, Oliveira AC de. Analysis of the air/water channels of gastrointestinal endoscopies as a risk factor for the transmission of microorganisms among patients. *Am J Infect Control*. 2012;40(10):913-916.
13. Ofstead CL, Wetzler HP, Doyle EM, et al. Persistent contamination on colonoscopes and gastroscopes detected by biologic cultures and rapid indicators despite reprocessing performed in accordance with guidelines. *Am J Infect Control*. 2015;43(8):794-801.
14. Cottarelli A, Giusti MD, Solimini AG, et al. Microbiological surveillance of endoscopes and implications for current reprocessing procedures adopted by an Italian teaching hospital. *Ann Di Igiene Medicina Prev E Di Comunita*. 2020;32(2):166-177.
15. Johani K, Hu H, Santos L, et al. Determination of bacterial species present in biofilm contaminating the channels of clinical endoscopes. *Infect Dis Heal*. 2018;23(4):189-196.
16. Kovaleva J, Meessen N, Peters F, et al. Is bacteriologic surveillance in endoscope reprocessing stringent enough? *Endoscopy*. 2009;41(10):913-916.
17. Ofstead CL, Smart AG, Hopkins KM, Wetzler HP. The utility of lighted magnification and borescopes for visual inspection of flexible endoscopes. *Am J Infect Control*. 2023;51(1):2-10.
18. Primo MGB, Tipple AFV, Costa D de M, et al. Biofilm accumulation in new flexible gastroscopes channels in clinical use. *Infect Control Hosp Epidemiology*. 2022;43(2):174-180.
19. Alfa MJ. Biofilms on instruments and environmental surfaces: Do they interfere with instrument reprocessing and surface disinfection? Review of the literature. *Am J Infect Control*. 2019;47:A39-A45.